


DESTINO™ TWIST

Deflectable Steerable Guiding Sheath

ORDERING GUIDE — Interventional Models


Product Description

The Oscor Destino™ Twist is a steerable guiding sheath that offers a deflectable tip with an ergonomic handle. The Destino™ Twist is the ultimate tool in gaining access to the most difficult to reach sites while maintaining hemostasis with Oscor SureSeal™ technology. The Destino™ Twist is your one sheath solution intended to facilitate the intracardiac, renal, and peripheral placements of diagnostic and therapeutic devices.



Features

- Ergonomic handle with clear imprint of the curve deflection, french size and sheath length
- Tip position indicator showing neutral or deflected curve
- Accurate tip deflection response to handle collar rotation
- Soft, atraumatic tip reduces potential for trauma and provides smooth transition to catheter
- Radiopaque tip marker enhances fluoro visibility
- Curve retention for accurate catheter placement
- Hydrophobic coated sheath improves access
- Braided Flexsteer™ shaft provides exceptional control torqueability, and kink resistance
- Snap-locking dilator connector
- Accepts maximum guidewire up to 0.038"/0.97 mm



P: 727.937.2511
F: 727.934.9835
E: sales@oscor.com
www.oscor.com

CAUTION: Federal Law restricts this device to sale by or on the order of a physician (U.S.A.). Please refer to the instructions for use for a complete listing of the indications, contraindications, precautions and warnings, where applicable.



DESTINO™

TWIST


Deflectable Steerable Guiding Sheath

ORDERING GUIDE — Interventional Models

Model # Set <i>(One per box)</i>	Sheath Inner Diameter	Usable Sheath Length	Usable Dilator Length	Curve Bend
STANDARD				
DST0654509	● 6.5 F	45 cm	63 cm	9 mm
DST0654517	● 6.5 F	45 cm	63 cm	17 mm
DST0654522	● 6.5 F	45 cm	63 cm	22 mm
DST0655507D2*	● 6.5 F	55 cm	73 cm	7 mm
DST0655509	● 6.5 F	55 cm	73 cm	9 mm
DST0655517	● 6.5 F	55 cm	73 cm	17 mm
DST0655522	● 6.5 F	55 cm	73 cm	22 mm
DST0656509	● 6.5 F	65 cm	83 cm	9 mm
DST0659009	● 6.5 F	90 cm	108 cm	9 mm
TD				
DST0654509TD	● 6.5 F	45 cm	63 cm	9 mm
DST0655509TD	● 6.5 F	55 cm	73 cm	9 mm
DST0656509TD	● 6.5 F	65 cm	83 cm	9 mm
STANDARD				
DST0704509	● 7 F	45 cm	63 cm	9 mm
DST0704517	● 7 F	45 cm	63 cm	17 mm
DST0704522	● 7 F	45 cm	63 cm	22 mm
DST0705509	● 7 F	55 cm	73 cm	9 mm
DST0705517	● 7 F	55 cm	73 cm	17 mm
DST0705522	● 7 F	55 cm	73 cm	22 mm
DST0706509	● 7 F	65 cm	83 cm	9 mm
DST0706530	● 7 F	65 cm	83 cm	30 mm
DST0709009	● 7 F	90 cm	108 cm	9 mm
TD				
DST0704509TD	● 7 F	45 cm	63 cm	9 mm
DST0705509TD	● 7 F	55 cm	73 cm	9 mm
DST0706509TD	● 7 F	65 cm	83 cm	9 mm

Model # Set <i>(One per box)</i>	Sheath Inner Diameter	Usable Sheath Length	Usable Dilator Length	Curve Bend
STANDARD				
DST0854517	● 8.5 F	45 cm	63 cm	17 mm
DST0855517	● 8.5 F	55 cm	73 cm	17 mm
DST0856517	● 8.5 F	65 cm	83 cm	17 mm
DST0856530	● 8.5 F	65 cm	83 cm	30 mm
DST0859017	● 8.5 F	90 cm	108 cm	17 mm
STANDARD				
DST1209030	● 12 F	90 cm	108 cm	30 mm
DST12010035	● 12 F	100 cm	118 cm	35 mm
STANDARD				
DST1357139	● 13.5 F	71 cm	89 cm	39 mm

Shaft Configurations



Set Contains

- (1) Destino™ Twist steerable guiding sheath
- (1) Destino™ Twist straight curved dilator

* Contains second soft dilator with 80 cm usable length.



P: 727.937.2511
F: 727.934.9835
E: sales@oscor.com
www.oscor.com

CAUTION: Federal Law restricts this device to sale by or on the order of a physician (U.S.A.). Please refer to the instructions for use for a complete listing of the indications, contraindications, precautions and warnings, where applicable.



DSTOGIN-112116-00

DESTINO™


TWIST

Deflectable Steerable Guiding Sheath

ORDERING GUIDE — Transseptal Models

Product Description

The Oscor Destino™ Twist is a steerable guiding sheath that offers a deflectable tip with an ergonomic handle. The Destino™ Twist is the ultimate tool in gaining access to the most difficult to reach sites while maintaining hemostasis with Oscor SureSeal™ technology. The Destino™ Twist is your one sheath solution intended to facilitate the intracardiac, renal, and peripheral placements of diagnostic and therapeutic devices.



Features

- Ergonomic handle with clear imprint of the curve deflection, french size and sheath length
- Tip position indicator showing neutral or deflected curve
- Accurate tip deflection response to handle collar rotation
- Soft, atraumatic tip reduces potential for trauma and provides smooth transition to catheter
- Radiopaque tip marker enhances fluoro visibility
- Curve retention for accurate catheter placement
- Hydrophobic coated sheath improves access
- Braided Flexsteer™ shaft provides exceptional control torqueability, and kink resistance
- Snap-locking dilator connector
- Accepts maximum guidewire up to 0.038"/0.97 mm



P: 727.937.2511
F: 727.934.9835
E: sales@oscor.com
www.oscor.com

CAUTION: Federal Law restricts this device to sale by or on the order of a physician (U.S.A.). Please refer to the instructions for use for a complete listing of the indications, contraindications, precautions and warnings, where applicable.




DESTINO™

TWIST

Deflectable Steerable Guiding Sheath

ORDERING GUIDE — Transseptal Models



Model # Set <i>(One per box)</i>	A1 Sheath Inner Diameter	A2 Usable Sheath Length	B Usable Dilator Length	A3 Curve Bend	C Compatible Transseptal Needle
8.5 F STANDARD					
DST0856717	● 8.5 F	67 cm	85 cm	17 mm	89 cm
DST0857317	● 8.5 F	73 cm	90 cm	17 mm	98 cm
DST0856722	● 8.5 F	67 cm	85 cm	22 mm	89 cm
DST0857322	● 8.5 F	73 cm	90 cm	22 mm	98 cm
DST0856750	● 8.5 F	67 cm	85 cm	50 mm	89 cm
DST0857350	● 8.5 F	73 cm	90 cm	50 mm	98 cm
12 F STANDARD					
DST1207122	● 12 F	71 cm	89 cm	22 mm	98 cm
DST1207522	● 12 F	75 cm	93 cm	22 mm	98 cm
DST1207139	● 12 F	71 cm	89 cm	39 mm	98 cm
DST1207150	● 12 F	71 cm	89 cm	50 mm	98 cm
DST1207550	● 12 F	75 cm	93 cm	50 mm	98 cm
13.5 F STANDARD					
DST1357139	● 13.5 F	71 cm	89 cm	39 mm	98 cm

Set Contains

- (1) Destino™ Twist steerable guiding sheath
- (1) Destino™ Twist straight dilator



P: 727.937.2511
F: 727.934.9835
E: sales@oscor.com
www.oscor.com

CAUTION: Federal Law restricts this device to sale by or on the order of a physician (U.S.A.). Please refer to the instructions for use for a complete listing of the indications, contraindications, precautions and warnings, where applicable.

